

PRODUCT SPECIFICATION

Product Name:	Hush-Slush™ Systems
Product Number:	ORS-1075EU Containing ORS-1075EU (Hush-Slush Machine) and pre-installed ORS-K24EU (Auto-Stir™). The Solution Warmer ORS-2038EU is part of the Hush-Slush System as well
Description:	Hush-Slush™ Machine, 230V, 7.5 Liters Round Basin
Qty/Cs:	1 unit per case
Commodity Code:	8543709099
Country of Origin:	US
Intended use:	The Hush-Slush System with its Auto-Stir™ technology is designed to cool surgical irrigation solutions and / or create surgical slush in order to cool delicate tissues. Next to that the Solution Warmer is designed to warm and maintain the temperature of surgical solutions prior to their use.

Specifications Hush-Slush Machine - ORS-1075EU

Power	220-240 Volt AC, 50/60 Hz, 3 Amps
Fuse	250 Volt AC; 5 Amp, 5 x 20 mm Time Delay (F5.0AL)
Temperature Range	-1 to 0°C
Med. Equip. Class.	Class I
EMC Compliance	The ORS-1075EU with ORS-2038EU and ORS-K24EU have been evaluated with respect to electric shock, fire and mechanical hazards in accordance with EN IEC 60601-1.
Equipment mass	50 Kg (including working load)
Dimensions	91 x 45 x 45 cm
Fluid volume	1000 ml (minimum) to 4000 ml (maximum)
Usage	Use only with Slush (and Warmer) Disc Drapes – Product Numbers ORS-130, ORS-320, ORS-321, ORS-325
Accessories Included	Auto-Stir™ – pre-installed, Product Number ORS-K24EU Non-sterile Wire Whisk, stainless steel, Product Number ORS-K12EU

Specifications Solution Warmer - ORS-2038EU

Power	220-240 Volt AC, 50/60 Hz, 0,75 Amps
Fuse	2x; 250 Volt AC; 1,25 Amp, 5 x 20mm Time Delay (F1.25AL)
Normal Temp Range	Room temperature to 49°C
Basin Temp Set Point	20 to 43 ± 1°C
Med. Equip. Class.	Class I
Equipment mass	6 Kg
Basin dimensions	29 x 15 x 10 cm
Fluid volume	1000 ml (minimum) to 3000 ml (maximum)
Usage	Use only with Slush (and Warmer) Disc Drapes – Product Numbers ORS-

	130, ORS-320, ORS-321, ORS-325
Included	Part of the Hush-Slush System

Specifications Auto-Stir - ORS-K24EU	
Power	120 Volt AC, 50/60 Hz, 0,6 Amps
Fuse	250 Volt AC; 0,8 Amp, 5 x 20mm Time Delay (F0.8AL)
Med. Equip. Class.	Class I
Included	Part of the Hush-Slush System

PRODUCT SPECIFICATION POWER CORDS

Product Name: Hush-Slush™
Product Number: SE-EU15
Description: Power Cord Kit 4.6 Meter EU
Qty/Cs: 1 unit per case
Dimensions: 23 x 23 x 10 cm
Weight: 700 gram
Packaging: Kit containing an EU Power Cord of 4.6 meter / 15ft

Product Name: Hush-Slush™
Product Number: SE-EU25
Description: Power Cord Kit 7.6 Meter EU
Qty/Cs: 1 unit per case
Dimensions: 23 x 23 x 10 cm
Weight: 900 gram
Packaging: Kit containing an EU Power Cord of 7.6 meter / 25ft

Product Name: Hush-Slush™
Product Number: SE-UK15
Description: Power Cord Kit 4.6 Meter UK
Qty/Cs: 1 unit per case
Dimensions: 23 x 23 x 10 cm
Weight: 700 gram
Packaging: Kit containing a UK Power Cord of 4.6 meter / 15ft

Product Name: Hush-Slush™
Product Number: SE-UK25
Description: Power Cord Kit 7.6 Meter UK
Qty/Cs: 1 unit per case
Dimensions: 23 x 23 x 10 cm
Weight: 900 gram
Packaging: Kit containing a UK Power Cord of 7.6 meter / 25ft

PRODUCT SPECIFICATION

Product Name:	Hush-Slush™ Systems
Product Number:	ORS-130, ORS-320, ORS-321, ORS-325
Description:	Disc Drapes Single use, sterile discs with integral drapes to fit into the Hush-Slush basin and the Solution Warmer
Components:	Polyurethane sheet, different sizes Polyurethane disc
Usage:	Single Use Use only with Hush-Slush™ Systems – Product Number ORS-1075EU
Qty/Cs:	24 units per case
Packaging:	Sterile
Sterilization method:	ETO

Part Number	Description	Dimensions	Box Dimensions	Box weight
ORS-130	Slush Disc-Drape	122 cm x 122 cm	52 x 33 x 25.4 cm	8 kg
ORS-320	Slush + Warmer Disc-Drape	168 cm x 112 cm	51 x 33 x 25.4 cm	10 kg
ORS-321	Slush + Warmer Disc-Drape	168 cm x 132 cm	51 x 33 x 25.4 cm	11.7 kg
ORS-325	Slush + Warmer Disc-Drape with Skirt	168 cm x 112 cm x 91 cm	56 x 45 x 94 cm	10 kg

Technical Properties	Performance
Package Seal Integrity	PASS
Puncture Testing	PASS
Biocompatibility	PASS
Color	Clear ORS-325 only: additional white non-woven skirt with elastic band